

# DINAMAP Pro 1000

Vital Signs Monitor



imagination at work





## The simple solution to accurate monitoring.

The DINAMAP® Pro 1000 monitor gives you core vital signs with 3-lead ECG and respirations in one simple yet comprehensive system. It is flexible for use throughout the hospital, from the ER to the GI lab, and is portable for transport wherever you need it.

The DINAMAP Pro 1000 has an easy-to-use design with an intuitive menu structure. Configurable features like spot SpO<sub>2</sub> enhance productivity and reduce time spent learning new technology.

Like all GE monitors, the DINAMAP Pro 1000 incorporates gold standard parameters to provide the most accurate and reliable vital signs readings. Its durable design ensures reliability so it is ready for use where you need it, when you need it.

And its easy connection to the Unity Network® lets you seamlessly maintain patient records throughout the enterprise, saving you time and allowing you to make faster, more appropriate care decisions.

# It works the way you do.

The DINAMAP Pro 1000 has an intuitive design that makes it easy to use throughout the hospital. Its configurable features include custom autocycles that automatically change the BP determination interval according to a user-defined program. And DINAMAP ASAP is not only fast, but adjusts inflation pressures for improved patient comfort and accurate results.


## NIBP Performance Improvements

SpO<sub>2</sub> spot check not only provides simple, fast, and accurate oxygen saturation, but also includes alarm technology that automatically returns the unit to “ready” status at the touch

of a button when doing routine spot checks. With your choice of SpO<sub>2</sub> technologies, the DINAMAP Pro 1000 enhances accuracy even in the most difficult patient situations, such as motion and low peripheral perfusion.

The DINAMAP Pro 1000 even helps prevent the loss of critical patient data by retaining information from one previous patient. And the standby key eliminates nuisance alarms and allows you to resume monitoring at a simple key press.

The large 10.4" color screen displays patient information on easy-to-see, customizable screens. Graphical and tabular trending speed retrieval of data, while the familiar Trim Knob® enhances menu access. Three simple default tables give you the information you need, when you need it so you can make the best patient care decisions.



Our latest DINAMAP ASAP algorithm reflects reductions in determination time and accuracy.



## Gold standards at market prices.

With the DINAMAP Pro 1000, you get unparalleled cardiac monitoring including gold standard parameters that provide reliable, accurate data so you can make the best care decisions faster:

- Legendary Marquette® ECG with EK-Pro™ and respirations offering comfort that lethal arrhythmias will be brought to your attention
- DINAMAP ASAP technology for accurate and reliable non-invasive BP monitoring
- Your choice of SpO<sub>2</sub> solutions – the industry leading Masimo® SET® (Gold Standard) or Nellcor OxiMax™
- Optional IVAC® Turbo
- Temp™ for temperature in seconds



## Easy connection to critical information.

The DINAMAP Pro 1000 quickly connects to the Unity Network with an easily-installed communications module. This connection allows caregivers to view DINAMAP Pro 1000 patient data on any network-connected workstation or device anywhere within or beyond the enterprise.

The compact communications module fully integrates into the back of the DINAMAP Pro 1000 monitor, providing one serial port and one Ethernet port for connecting to the Unity

Network. This simple device offers a cost effective solution for accessing DINAMAP Pro 1000 patient data at the central station, via the StatView™ alarm notification system, and on other GE monitoring products across the enterprise.

The Unity Network supports both wired and wireless viewing devices to ensure that life-critical data such as waveforms and monitor alarms are consistently, safely and securely delivered to caregivers wherever they are.



For more than 100 years, scientists and industry leaders have relied on General Electric for technology, services and productivity solutions.

So no matter what challenges your healthcare system faces – you can always count on GE to help deliver the highest quality services and support.

For details, please contact your GE representative today.

**GE Healthcare**

Waukesha, Wisconsin U.S.A.  
[www.gehealthcare.com](http://www.gehealthcare.com)

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

©2004 General Electric Company  
GE Medical Systems, a General Electric company, going to market as GE Healthcare.

04-8997 8/04 Printed in U.S.A.

imagination at work 